

A Benchmark Comparison

Novel membrane stirrer prevents excessive foaming in rhamnolipid production using *Pseudomonas putida*

Introduction

Biosurfactants as Rhamnolipids (RLs) are naturally occurring amphipathic surface-active agents produced by microorganisms as yeast and bacteria. RLs are eco-friendly alternatives to petrochemical-derived synthetic surfactants with wide-ranging applications across various industries including environment, petroleum, agriculture, food and pharmaceuticals.

Apart from the typical microbial cultivation challenges, producing biosurfactants is especially prone to extensive foaming. When gas bubbles are introduced into the fermentation broth, the biosurfactants adsorb at the gas-liquid interfaces, forming bubbles with micelles stabilizing the foam (see Fig. 1) [1]. Foaming in bioprocesses can reduce product quality and quantity (e.g., by protein denaturation) as well as productivity and result in a loss of the whole-cell biocatalyst [2-4]. In this application note, fermentation experiments using *Pseudomonas putida* KT2440 SK4 strains were conducted to establish a benchmark run comparing RL titers in the BioThrust MemStir with a conventional stirred-tank reactor (STR) at a 2 L scale. The MemStir demonstrated a space-time yield (STY) of up to $118 \text{ mg}_{\text{RL}} \text{ L}^{-1} \text{ h}^{-1}$, ranking it high in industrial RL production with both conventional and complex fermentation methods.

Materials & Methods




Figure 1: Overview of RL production with *P. putida* by either conventional sparged bubble-aeration in an STR (A) or membrane aeration with BioThrust's MemStir (B). Created with BioRender.com.

Objective

- Investigate foaming, oxygen supply solutions & transmembrane pressure (TMP) in RL production and its impact on product quantity.
- Performance analysis of BioThrust's MemStir with a conventional STR in RL batch processes using *Pseudomonas putida*.
- Getting an insight into a complex RL production using a fed-batch process.

Key Results

1. Foam-free batch and fed-batch state-of-the-art fermentation.
2. More than 3-fold increase of rhamnolipid production.
3. Increased productivity due to 65 % increased batch volume.
4. The MemStir enabled a superior STY of up to $118 \text{ mg}_{\text{RL}} \text{ L}^{-1} \text{ h}^{-1}$.

Specifications

- **Bacterial strains:**
P. putida KT2440 SK4 and *P. putida* KT2440 Δ flag SK4
- **Culture medium:**
 - Preculture 1: Lysogeny Broth (LB).
 - Preculture 2: 3x buffered Minimal Salt Medium (MSM) with 10 g L^{-1} glucose.
- **Bioreactor:**
Working volume of 1.2 L (STR) and 2 L (MemStir).
- **Process modalities:**
Batch and Fed-Batch
(DO controlled feed: 200 g L^{-1} glucose).
- **Process time:**
Batch: 10 h, Fed-Batch: 18 h.

Abbreviations

RL = Rhamnolipid	CDW = Cell Dry Weight
STR = Stirred-tank reactor	XO ₂ = O ₂ content in the feed gas
STY = Space-time yield	MSM = Minimal salt medium
SK4 = Wild type containing RL synthesis promoters	DO = Dissolved oxygen
LB = Lysogeny broth	GLC = Glucose
	rpm = Revolutions per minute
	TMP = Transmembrane pressure

Parameter	MemStir	STR
Stirrer	100 – 350 rpm	300 – 1200 rpm (Rushton turbines)
DO	>30 %	>30 %
pH	7.0 (Base adjusted)	7.0 (Base adjusted)
Temperature	30 °C	30 °C
Aeration	1 L/min	0.12 L/min
TMP	0 – 0.3 bar	–
Antifoam	–	23 g/L

Results

For the first time a novel membrane impeller, the MemStir was developed and patented (WO 2021/152128) and applied for a benchmark run for a foam-free RL production (see Fig. 1). Please refer to the Technical Application Note for further details. The benchmark fermentation was conducted in a conventional STR using submerged bubble-aeration and a high titer RL-producing recombinant *P. putida* KT2440 Δ flag strain [5]. The MemStir batch fermentation was conducted in biological duplicates with *P. putida* KT2440 SK4 for 10 h. The MemStir replaced the stirrer and ring sparger as the mixing and aeration unit for bubble-free aeration of the bioreactor. Figure 3A and 3C show the cell dry weight (CDW), glucose (Glc), ammonium sulfate ((NH₄)₂SO₄), and RL concentrations within the conventional bioreactor process (A) and within the MemStir process (C). Figures 3B and D show the DO and O₂ content in the feed gas X_{O2} as well as the stirring rate within the conventional bioreactor (B) and the MemStir (D).

Batch Benchmark




Figure 2: Images of batch fermentations of *P. putida* in BioThrust's MemStir (A) and in an STR with conventional stirring and aeration (B). Working volumes of the conventional STR experiments were reduced due to occurring foam formation.




Figure 3: Results of the batch cultivations with *P. putida* KT2440 Δ flag SK4 in the conventional STR (A, C) and with *P. putida* KT2440 SK4 in the MemStir STR (B, D). (A) and (C) are displaying the cell dry weight (CDW), glucose (Glc), ammonium sulfate ((NH₄)₂SO₄), and rhamnolipid (RL) titers over time. (B) and (D) indicate the DO, stirring rate and O₂ content in the feed gas X_{O2} over time. Experiments were conducted as biological duplicates, whereas off-line samples were analyzed as technical triplicates.

In the MemStir a maximum biomass concentration of 2.57 g L^{-1} with a $\mu_{\max}$ of 0.59 h^{-1} was achieved. The DO was maintained at or above 30 % throughout the cultivation. At the end of the fermentation, a RL titer of 1.17 g L^{-1} was achieved, which is a measured 3.5-fold higher than in the benchmark STR fermentation using bubble-aeration. The culture volume was 1.7-fold higher in contrast to the benchmark fermentation. Thus, higher product amounts per batch can be achieved using the MemStir. Neither foam formation nor oxygen limitation were observed over the fermentation time. More RL accumulated in the MemStir process, increasing the titer approximately 3.3-fold with an STY of $118 \text{ mg}_{\text{RL}} \text{ L}^{-1} \text{ h}^{-1}$.

Fed-Batch investigation




Figure 4: Fed-Batch cultivation of *P. putida* KT2440 SK4 in the MemStir setup ($n = 2$) displaying CDW, Glc, (NH₄)₂SO₄ and RL titers (E) as well as DO, stirring rate and X_{O2} (F). Off-line samples were analyzed as technical triplicates.

The DO-based fed-batch fermentation was conducted in biological duplicates with *P. putida* KT2440 SK4 (see Fig. 4). A concentrated MSM containing 200 g L^{-1} glucose was used for feeding. Fermentation began with a batch phase that concluded after glucose depletion at 9 hours. Following depletion, DO levels rose to 70 %, prompting the start of the feed phase (Fig. 4F). At the end of the fermentation, 5.8 g rhamnolipids (2.6 g L^{-1}) and 14.1 g CDW (6.3 g L^{-1}) were reached (Fig. 4E). Compared with the batch fermentation with the SK4 strain, the RL titer (2.6 g L^{-1}) was increased by 2.2-fold. Furthermore, a STY of $118 \text{ mg}_{\text{RL}} \text{ L}^{-1} \text{ h}^{-1}$ was achieved, equaling to the batch fermentation.

Discussion

The benchmark experiment comparing BioThrusts MemStir-driven bioreactor with a conventional STR demonstrated several advantages:

1. Due to no foam formation, the working volume could be increased to the maximum value recommended by the reactor manufacturer, which in turn allows for a more accurate assessment of the reactor's productivity.
2. RL bioprocesses will no longer require a continuous supply of antifoam, reducing downstream processing costs, purification expenses, and potential product loss [6].
3. BioThrust's MemStir achieved a more than 3-fold higher rhamnolipid titer compared to the conventional STR. In the MemStir aerated batch as well as fed-batch fermentations, a STY up to $118 \text{ mg}_{\text{RL}} \text{ L}^{-1} \text{ h}^{-1}$ was reached.

The MemStir bioreactor demonstrates the feasibility of prolonged operation with minimal risk of batch failure, thereby increasing the potential yield. Both batch and fed-batch processes exhibited product titers and STY that are competitive with those reported in recent literature.

Conclusion

In this study, also published by Bongartz et al., (2023), high-titer RLs were produced by bubble-free aeration, preventing any foam formation and oxygen limitation during the fermentation process. Compared to previous reports about recombinant rhamnolipid batch production with *P. putida* KT2440, obtained product titers were greatly increased to conventional processes. Overall, cell concentrations and STY can be ranked in a higher performance category. While there are other process strategies that increase the system effectivity even further, results of the fed-batch shows a promising potential to further exceed current industrial standards when optimizing operation and feeding strategies.

References

- [1] – Junker (2007). *Foam and its mitigation in fermentation systems*. *Biotechnology Progress*. 23(4), pp. 767–784. DOI: 10.1021/bp070032r.
- [2] – Vardar-Sukan (1998). *Foaming: Consequences, prevention and destruction*. *Biotechnology Advances*. 16(5–6), pp. 913–948. DOI: 10.1016/S0734-9750(98)00010-X.
- [3] – Demling et al. (2020). *Selection of a recyclable in situ liquid–liquid extraction solvent for foam-free synthesis of rhamnolipids in a two-phase fermentation*. *Green Chemistry*. 22(23), pp. 8495–8510. DOI: 10.1039/D0GC02885A.
- [4] – Handa-Corrigan et al. (1989). *Effect of gas-liquid interfaces on the growth of suspended mammalian cells: Mechanisms of cell damage by bubbles*. *Enzyme and Microbial Technology*. 11(4), pp. 230–235. DOI: doi.org/10.1016/0141-0229(89)90097-5.
- [5] – Tiso, T. et al. (2020) 'Integration of Genetic and Process Engineering for Optimized Rhamnolipid Production Using *Pseudomonas putida*', *Frontiers in Bioengineering and Biotechnology*, 8, p. 565499. Available at: <https://doi.org/10.3389/FBIOE.2020.00976/BIBTEX>.
- [6] – Anic et al. (2018). *Production of rhamnolipids by integrated foam adsorption in a bioreactor system*. *AMB Express*. 8(1), p. 122. DOI: 10.1186/s13568-018-0651-y.
- [7] – Bongartz et al. (2023). *A novel membrane stirrer system enables foam-free biosurfactant production*. *Biotechnology & Bioengineering*. 120(5), p. 1282, Table 1. DOI: 10.1002/cite.202255100.



BioThrust GmbH
Pauwelsstraße 17
52074 Aachen
Germany

Dr.-Ing. Patrick Bongartz
patrick.bongartz@biothrust.de

Authors: Tobias Karmainski, Dr. Yasemin van Heuvel, Tim Milfeit, Moritz Meyer,
Prof. Dr.-Ing. Matthias Wessling, Dr. Till Tiso, Prof. Dr. Lars Mathias Blank,
Dr.-Ing. Patrick Bongartz

Copyright ©2024 BioThrust GmbH. All rights reserved.



www.biothrust.com